

ComforTrac

comfortable • simple • effective

Lumbar Traction



The ComforTrac lumbar home traction unit was carefully engineered to provide maximum comfort, efficacy, and ease of use. Details that highlight this superior design include an automatic carriage return and innovative hinged design that requires no set-up. Likewise, transporting the unit has been made effortless through the use of lightweight durable materials and a carry-case with smooth rolling in-line skate wheels and a retractable stow-away handle typically found only in top-quality ballistic luggage.

Lumbar Traction

Our innovative design makes it possible for patients to go home with the same quality traction they receive during physical therapy.

Memory Pillow
conforms comfortably to the head

Padded Surface
provides superior comfort

Unit folds
in half for easy storage and transportation



Hand Pump
our patented hand pump makes it easy to increase or decrease traction force

Quick Release Button
push to release the traction force

Smooth Traction
provides consistent and even force

SPECIFICATIONS:

Maximum force: 200 lbs.

Hand pump: Displays force in both pounds and kilograms

Headrest material: Memory foam

Traction bed material: Cushion foam

Testing: 5,000 cycles which equals using 3 times daily for 5 years

SUPPORTING RESEARCH AND REFERENCES:

Komori H, Shinomiya K, Nakai O, et. al., The Natural History of Herniated Nucleus Pulposus with Radiculopathy. Spine 21:225-229, 1996.

Onel D, Tuzlaci M, Sari H, et. al. Computed Tomographic Investigation of the Effect of Traction on Lumbar Disc Herniations. Spine 14:82-90, 1989.

Matthews JA, Hickling J, Lumbar Traction: A Double-Blind Controlled Study for Sciatica. Rheumatology and Rehabilitation 14:222-225, 1975.



ComforTrac Inc.
3930 Horseshoe Bend Road
Keyville, VA 23947

866.378.4672
703.891.0455
703.288.3304 fax

info@comfortrac.net
www.comfortrac.net